

## RVIT SERIES



Elements of **RVIT series carts** are made of:

**frame** of aluminum profile, powder coated in colour consistent with TECH-MED RAL palette

**filling of the cabinet** (optional): powder coated steel in white, stainless steel grade 0H18N9 (AISI 304), aluminum or ABS plastic in white

**top** (optional): powder coated steel in white, stainless steel grade 0H18N9 (AISI 304), aluminum or ABS plastic in white, with bands on three side

**base** (optional): powder coated steel in white, steel with white ABS plastic cover or stainless steel grade 0H18N9 (AISI 304), equipped with 125 mm diameter castors, including two with lock

**decorative elements** powder coated in colour consistent with TECH-MED RAL palette

**Option:** ABS version - there is a possibility of making above white elements in gray

**Equipment and additional accessories:** pull-out side top, tilting containers, shelf, basket for accessories, waste bowl, instrument container, glove box container, used needles container, catheter basket, waste bin, IV drip pole, power strip, lock for drawer, central lock

**Dimensions of cabinet:**

length: **600 mm, 700 mm, 800 mm** [L]

width: **500 mm** [S]

height from the ground to the top: **1000 mm** [H]

height of cart with top section: **1700 mm** [H]

RVIT-511-B(ST)S(ST)F(ST)P(ST)-3000-K13-W

Examples of carts:



RVIT-1D11-B(ABS3)S(ST)F(ST)P(ABS3)-3000-K12-W



RVIT-2W01-B(KO)S(KO)F(KO)P(KO)-3000-K13-W



RVIT-301-B(ABS3)S(ST)F(ST)P(ST)-6018-K13-W



RVIT-401-B(ST)S(ST)F(ST)P(ST)-1018-K13-W



RVIT-401-B(ABS6)S(ST)F(ST)P(ABS6)-5012-K13-W



RVIT-501-B(ABS6)S(ST)F(ST)P(ABS6)-9005-K13-W