



HECO-01



HSM-01

Disinfection Stands

HECO-01

Disinfection stand HECO-01 is made of:

frame: powder-varnished aluminum profile in white with possibility for installation of additional equipment

base: powder-varnished steel equipped with feet with possibility of levelling

holder for disinfectant liquid: stainless steel, adjustable in height in a range of 90-210 mm for disinfectant liquid container (height of bottle without dispenser)

Additional equipment: gloves container, accessories basket, waste bag holder, waste bin with swing lid, information board

Dimensions: 500x460x1700 mm

HSM-01

Disinfection stand HSM-01 is made of:

base: plastic, equipped with feet with possibility of levelling

column: stainless steel

holder for disinfectant liquid: stainless steel, adjustable in height in a range of 90-210 mm for disinfectant liquid container (height of bottle without dispenser)

Additional equipment: container for gloves, accessories basket

Base diameter: 575mm

Height: 1300mm

* It is possible to modify the product according to the individual customers' needs.



HECO-01